

AquaTuff™ Wrap&Stow™ Thermocouple Instruments

Wrap&Stow™ designs are available with a unique, cable storage channel so that the heavy duty, patented probe can be stored safely alongside the unit housing. The Wrap&Stow™ insertion probe is factory-calibrated with the instrument for a higher degree of total system accuracy.

	35132	35135	35140	35232	35235	35240	35340
Probe Included:	DuraNeedle	Bell Surface	MicroNeedle	DuraNeedle	Bell Surface	MicroNeedle	MicroNeedle
Temperature Range:	-100° to 500°F (-73° to 260°C)	-100° to 500°F (-73° to 260°C)	-100° to 500°F (-73° to 260°C)	-100° to 500°F (-73° to 260°C)	-100° to 500°F (-73° to 260°C)	-100° to 500°F (-73° to 260°C)	-100° to 500°F (-73° to 260°C)
Accuracy:	±0.9°F (±0.5°C) total system accuracy	±0.5°F (±0.3°C)*	±0.9°F (±0.5°C) total system accuracy	±0.9°F (±0.5°C) total system accuracy	±0.5°F (±0.3°C)*	±0.9°F (±0.5°C) total system accuracy	±0.9°F (±0.5°C) total system accuracy
Resolution:	0.1°	0.1°	0.1°	0.1°/1° selectable	0.1°/1° selectable	0.1°/1° selectable	0.1°
Housing:	ABS Plastic	ABS Plastic	ABS Plastic	ABS Plastic	ABS Plastic	ABS Plastic	ABS Plastic
Hold:	No	No	No	Yes	Yes	Yes	-
Backlight:	No	No	No	Yes	Yes	Yes	-
Power:	(2) 1.5V AAA	(2) 1.5V AAA	(2) 1.5V AAA	(2) 1.5V AAA	(2) 1.5V AAA	(2) 1.5V AAA	(2) 1.5V AAA
Battery Life:	1800 hours	1800 hours	1800 hours	1800 hours	1800 hours	1800 hours	1800 hours
Auto Off:	10 min.	10 min.	10 min.	10 min.	10 min.	10 min.	10 min.
Weight:	7 oz (199 g)	8 oz (227 g)	7 oz (199 g)	7 oz (199 g)	8 oz (227 g)	7 oz (199 g)	7 oz (199 g)
Regulatory Listings:	CE (NSF) RoHS	CE RoHS	CE (NSF) RoHS	CE (NSF) RoHS	CE RoHS	CE (NSF) RoHS	CE RoHS
Warranty:	5 Year	5 Year	5 Year	5 Year	5 Year	5 Year	5 Year

* Accuracy specification for instrument only.

35132

AquaTuff™ Wrap&Stow™
Thermocouple Instrument
with DuraNeedle Probe



35235

AquaTuff™ Wrap&Stow™
Thermocouple Instrument
with Surface Probe



35340

AquaTuff™ Wrap&Stow™ ITS™
Thermocouple Instrument
with MicroNeedle Probe



Instrument with ITS™ and Memory Storage



The Intelligent Temperature Stabilization (ITS™) feature, on the 35340 AquaTuff™ Wrap&Stow™ ITS™ Thermocouple, prevents the temperature from being displayed until the final stabilized temperature is reached. While in the ITS mode, you have the option of recording the stabilized temperature into the 35340 memory. The memory can store up to 250 readings, which can then be reviewed by scrolling up or down.